



There is no annual preventative maintenance required for Hopkins Medical Products Pulse Oximeters.

**Instructions for Cleaning and Disinfecting Hopkins Medical Products Pulse Oximeters:**

- Surface-clean sensor with a soft gauze by wetting with a solution such as 70% isopropyl alcohol, if low-level disinfection is required, use a 1:10 bleach solution.
- Then Surface-clean with clean water and wipe down with a dry cloth.

**Caution:** Do not sterilize by steam, or ethylene oxide. Do not use the sensor if it is damaged.

2017 Hopkins Medical Products